



## Justice Mark II

### Rugged Protective Case

#### Durable Protection for Your NONIN Finger Pulse Oximeter

The Justice Mark II is a lightweight, rugged case designed to safeguard the Nonin finger pulse oximeter in extreme environments. The case has an automatic purge valve to compensate for dramatic changes in pressure and is watertight for maximum protection in almost any situation.



The case offers a crushproof exterior to keep the finger pulse oximeter safe and clean while out in the field. The customized interior easily stores the oximeter, lanyard, and extra batteries for your convenience.

Nonin developed the Justice Mark II to meet the needs of the medics of the U.S. Special Forces, working in the most punishing environments. With the Justice Mark II you can rest assured your Nonin finger pulse oximeter will be protected wherever you go.

### Product Highlights

- ✓ Durable - rugged crushproof design
- ✓ Watertight - leak proof to 50 ft./15 meters
- ✓ Automatic Purge Valve - pressure compensation
- ✓ Compartmentalized - easily stores a finger oximeter, lanyard, and extra batteries

### Ordering Information

|                      | Model # | Part #   |
|----------------------|---------|----------|
| Justice Mark II Case | JM2-CC  | 3802-000 |



[www.nonin.com](http://www.nonin.com)

Nonin Medical, Inc.  
 13700 1st Avenue North  
 Plymouth, MN 55441-5443, USA  
 Phone: 763.553.9968 Toll Free: 800.356.8874  
 Fax: 763.553.7807  
 info@nonin.com

✓ Where People Make the Difference